CareSens[®] H Beat [®] Blood Glucose Monitoring System

User Manual

Table of Contents

01	Important Information: Read This First!2
02	Specifications3
03	CareSens H Beat BT Blood Glucose Monitoring System4
04	Inserting or Replacing the Battery5
05	Caring for Your System6
06	CareSens H Blood Glucose Test Strip6
07	CareSens H Beat BT Blood Glucose Meter8
08	CareSens H Beat BT Blood Glucose Meter Display9
09	Setting Up Your System10
10	Pairing your Meter with Smartphone15
11	Setting the Strip Expiration Date Indicator17
12	Checking the System17
13	Using the Lancing Device20
14	HI and Lo Messages26
15	Target Blood Glucose Ranges26
16	Transferring Test Results Using Cable27
17	Meter Memory27
18	Setting the Alarm Function29
19	Understanding Error Messages31
20	General Troubleshooting33
21	Performance Characteristics34
22	Warranty Information36

PGUAA-0000304 REV1 2023-02

01 Important Information: Read This First!

For optimum safety and benefits, please read the entire manual contents before using the system.

Intended use:

The CareSens H Beat BT Blood Glucose Monitoring System is intended for the quantitative measurement of alucose in fresh capillary whole blood from the fingertip. The CareSens H Beat BT Blood Glucose Monitoring System is intended for use outside the body (in vitro diagnostic use) and is intended for use as an aid to monitor the effectiveness of diabetes management. The system is for self-testing or healthcare professional use and should not be used for the diagnosis of or screening for diabetes. In clinic and hospital settings, venous, arterial, and neonatal whole blood may also be used to measure blood glucose when drawn by healthcare professionals

Meaning of Symbols Used:

European Union

- CE Mark Authorised representative in
- Consult instructions for use WEEE (waste electrical and R **EC REP** the European Community/ electronic equipment) Temperature limit Manufacturer

Serial number

Use-by date

- IVD In vitro diagnostic medical device Caution LOT Batch code
- 🗞 Biological risks
- On the second second
- Glucose in blood samples reacts with the chemical in the test strip to produce a small electrical current. The CareSens H Beat BT meter detects this electrical current and measures the amount of glucose in the blood sample.
- The CareSens H Beat BT blood glucose meter is designed to minimise code related errors in monitoring by using the no-coding function
- The CareSens H Beat BT blood glucose meter should be used only with the CareSens H test strips.
- An abnormally high or low red blood cell count (hematocrit level
- over 70 % or below 15 %) may produce inaccurate results. • If your test result is below 60 mg/dL (3.3 mmol/L) or above 240 mg/dL (13.3 mmol/L), consult a healthcare professional
- immediately. · Inaccurate results may occur in severely hypotensive individuals or
- patients in shock. Inaccurate low results may occur for individuals experiencing a hyperglycemic- hyperosmolar state, with or without ketosis. Critically ill patients should not be tested with blood alucose meters.

• Do not use during or within 24 hours of receiving xylose absorption testing as it may cause inaccurate results.

If you need assistance, please contact your authorised i-SENS sales representative or visit www.i-sens.com for more information.

02 Specifications

Product specifications	;
Measurement range	10–600 mg/dL (0.6–33.3 mmol/L)
Sample size	Minimum 0.5 µL
Test time	5 seconds
Sample type	 Fresh capillary whole blood Healthcare professional use only: Fresh venous whole blood Fresh arterial whole blood Fresh neonatal whole blood
Calibration	Plasma-equivalent
Assay method	Electrochemical
Battery life	1,000 tests
Power	One 3.0 V lithium battery (disposable, type CR2032)
Memory	1,000 test results
Size	95 x 49.6 x 17.9 mm
Weight	52.8 g (with battery)

Bluetooth technology	 Frequency range: 2.4–2.4835 GHz Operating range distance: maximum 10 meters (unobstructed) Operating channels: 40 channels Security encryption: 128-bit AES (Advanced Encryption Standard) 			
Operating ranges				
Temperature		5–45 °C		
Relative humidity		10–90 %		
Hematocrit		15–70 %		
Storage/transport Co	nditio	ns		
		Meter (with battery)	0–50 °C	
Temperature		Test strip	1–30 °C	
		Control solution 8–30 °C		
Relative humidity		Test strip	10–90 %	

Relative humidity Test strip 03 CareSens H Beat BT Blood Glucose Monitoring

System

CareSens H Beat BT Blood Glucose Monitoring System includes the following items:

- * CareSens H Beat BT Blood Glucose Meter
- * User Manual
- CareSens H Beat BT Blood Glucose Monitoring System may include following items
- * CareSens H Blood Glucose Test Strips * Lancets * Lancing Device * Quick Reference Guide ' Battery * Loabook
- Carrying Case Check all the components after opening the CareSens H Beat BT Blood Glucose Monitoring System package. The exact contents are
- listed on the main box. • The cable for data management software can be ordered separately. Please contact your authorised i-SENS sales representative. 4
- 04 Inserting or Replacing the Battery
- The CareSens H Beat BT meter uses one 3.0 V lithium battery. Before using the meter, check the battery compartment and insert a battery if empty.

When the symbol appears on the display while the meter is in use, the battery should be replaced as soon as possible. The test results may not be saved if the battery runs out.



Make sure the meter is turned off. Push the cover in the direction of the arrow to open the battery compartment.



Step 2 Remove the used battery. Slip your index finger under the battery to lift and pull out as shown. Insert one new battery with the + side facing up and make sure the battery is inserted firmly.

Place the cover on the battery compartment. Push it down until you hear the tab click into place.

O Note

2

Removing the meter battery will not affect your stored results. However you may need to reset your meter settings. See page 10. 5

05 Caring for Your System

- Use a soft cloth or tissue to wipe the meter exterior. If necessary, dip the soft cloth or tissue in a small amount of alcohol. Do not use organic solvents such as benzene or acetone, or household and industrial cleaners that may cause irreparable damage to the meter. Caution:
- Do not expose the meter to direct sunlight, heat, or excessive humidity for an extended period of time. It is recommended to store and use the test system indoors.
- Do not let dirt, dust, blood, or water enter into the meter's test strip • Do not drop the meter or submit it to strong shock.
- Do not try to fix or alter the meter in any way.
- Avoid getting any liquid or moisture in the test strip vial. This can
- affect the test strips and cause inaccurate test results • Store all the meter components in the carrying case to prevent loss and help keep the meter clean.
- Do not apply samples other than capillary, venous, arterial, neonatal whole blood or control solution to the test strip.
- Store the meter in a cool and dry place between 0–50 °C.

Disposal of the meter If you need to throw your meter away, you should follow existing policies and procedures of your own country or region. For information about correct disposal, please contact your local council or authority. If you need assistance, contact your authorised i-SENS sales representative or visit <u>www.i-sens.com</u>.

06 CareSens H Blood Glucose Test Strip

The CareSens H Beat BT Blood Glucose Monitoring System measures blood glucose quickly and accurately. It automatically absorbs the small blood sample applied to the edge of the strip.


\bigcirc	
I.I.	 Contact bars Gently push the test strip, with its contact bars facing up, into the test strip port of meter
CareSe	Confirmation window Check here to see whether sufficient blood sample has been applied
	Edge to apply blood sample Apply blood sample here for testing

Warning!

- The CareSens H test strips should be used with fresh capillary whole blood samples, or with fresh venous, neonatal and arterial whole blood samples if drawn by healthcare professionals. Besides whole blood samples, serum or plasma samples can affect test results.
- · Fresh venous and arterial whole blood specimens containing the anticoagulants EDTA and Heparin are acceptable. Iodoacetate or fluoride/oxalate should not be used.
- Neonatal capillary samples may be drawn from the heel stick, not from the neonatal cord blood samples.
- Do not reuse test strins.
- Do not use test strips past the expiration date.
- Test strips in new, unopened vials and test strips in vials that have been opened can be used up until the expiration date printed on the test strip box and vial label if the test strips are used
- and stored according to its storage and handling methods. Store test strips in a cool and dry place at a temperature between 1–30 $^\circ$ C and 10–90 %relative humidity.
- Keep test strips away from direct sunlight or heat and do not freeze.
- Store test strips only in their original vial.
 Close the vial tightly after taking out a test strip for testing and use the strip immediately. Avoid getting any liquid or moisture in the test strip vial. This can affect the test strips and
- cause inaccurate test results. Do not apply samples other than capillary, venous, arterial, neonatal whole blood or control
- solution to the test strip. Handle test strips only with clean and dry hands.
- Do not bend, cut, or alter test strips in any way.
- For detailed storage and usage information, refer to the CareSens H test strip package insert.

▲ Caution

• Keep the meter and testing supplies away from young children. • Drying agents in the vial cap may be harmful if inhaled or swallowed and may cause skin or eye irritation.



08 CareSens H Beat BT Blood Glucose Meter Display


(1) **Memory symbol:** appears when test results stored in the memory are displayed] 🗘 🖡 🖇 🕬 🕯 (2) Setting symbol: appears when in SET mode


(3) PP2 alert: appears when the post-meal alarm has been set (4) Bluetooth: appears when Bluetooth feature is

Mute symbol: appears only when the sound

- is set to OFF 6 Control Solution flag: appears when the control
- solution test results are saved or displayed 7 Test results: test results displaying panel
- (8) KET symbol: appears when the test result is greater than 240 mg/dL
- (9) Hi: appears when the test result is greater than the selected hyperglycemia level
- (1) Smile symbol: appears when the test result is within the selected normal blood glucose range
- (1) Lo: appears when the test result lower than the selected hypoglycemia level (12) alarm: appears when the time alarm has been set
- mmol/L. ma/dL: unit for measuring blood glucose
- (14) Battery symbol: indicates meter battery is running low and needs to
- be replaced (15) Blood insertion symbol: indicates meter is ready for the application of
- a drop of blood or control solution (f) Fasting test flag: used for tests done after fasting for at least 8 hours
- Post-meal test flag: used for tests done after eating
- Pre-meal test flag: used for tests done before eating
- Month/Day/Hour/Minute: appears date and time
- O Note

t is recommended to check if the display screen on the meter matches the illustration above every time the meter turns on. Do not use the meter if the display screen does not exactly match the illustration as the meter may show incorrect results.





O Note

Step 10

screen displays 'OFF'.

09 Setting Up Your System

Press and hold the S button for 3 seconds to enter SET mode. After all settings are finished, press and hold the **S** button for 3 seconds to turn off the meter. Press the ◀ or ► button to change values. Press and hold the ◀ or ► button to scroll faster



ŚE

O Note

0 *

Setting Up Bluetooth

0 (1)

20/23

0

11

05-04 I2:00~

05-04 ¹2 h

05:04 12:00

0

ò

);;;00 i ib;;00~

 \leftrightarrow

Step 1 Entering the SET Mode Press and hold the S button for 3 seconds to enter SET mode. After all the segments flash across the screen, 'SET' will show up. Press the ◀ or ► button to select 'YES' and press the **S** button to go to the next step.

To turn on the Bluetooth feature.

The meter screen will display

press the ◀ or ▶ button.

'OFF' and 'On' in turn.

Follow steps 2–3 to pair your meter and smartphone. Pairing allows the meter to communicate wirelessly with your smartphone. Ensure that devices are within the naximum Bluetooth range (10 meters). Before pairing your meter and smartphone, download and install the SmartLog mobile app on your smartphone.



② Press the S button when 'On' blinks on the

screen. ③ If you do not want to connect your meter to your smartphone, press the S button when 'OFF' blinks on the screen

The meter will directly go to Step 3 Setting the Year on page 11. 10

Adjusting the Date and Time

/ Step 3 Setting the Year Press the ◀ or ► button to adjust until the correct year appears. When the present year appears, press the S button to confirm your selection and to go to the next

Step 4 Setting the Month

A number indicating the month will blink on the screen. Press the ◀ or ▶ button until the correct month appears. Press the S button to confirm your selection and to go to the next step.

Step 5 Setting the Date

Press the ◀ or ► button until the screen displays the correct date. Press the S button to confirm the date and to go to the next step.

Step 6 Setting the Time Format

The meter can be set in the AM/PM 12-hour or the 24-hour clock format.

Press the \triangleleft or \blacktriangleright button to select a format. The AM•PM symbol is not displayed in the 24-hour format. After selecting the format, press the S button to go to the next step.

Step 7 Setting the Hour

Press the < or > button until the correct hour appears. After the hour is set, press the S button to go to the next step.



Step 8 Setting the Minute

Press the ◀ or ► button until the correct minute appears. After setting the minute, press the S button to go to the next step.

On pressing the ◄ or ► button, the screen will display 'On' or 'OFF'. Press the S button to confirm the

> The meter will beep in the following instances if set to On. • When you push a button to turn on the meter, When the test strip is inserted in the meter.

• When the blood sample is absorbed into the test strip and the test starts.

 When the test result is displayed, • When you press and hold the < button to set the post-meal (PP2) alert


• When it is time for a pre-set blood glucose test. If the sound is set to OFF, none of the sound functions will work.

After setting the sound, press the **S** button to go to the next step.

The \mathbb{Z} symbol is displayed only when the sound is set to OFF. Turning on the Strip Expiration Date Indicator

This setting allows you to turn the strip expiration date indicator on or off. This setting turns the function on or off only. See page 17 to set the strip expiration date.

When 'EP' appears on the screen, press the ◀ or ► button. The screen will display 'On' or 'OFF'. Press the S button to confirm the setting. If you do not want to set the indicator, press the S button while the






when the test result is displayed right after the test. If the expiration date is set to October of 2025, the meter will display 'EP' at the start of November, 2025.

1101,42:37~

ILL



than the selected level. Press the ◀ or ► button until the desired hypoglycemia level between 10 and 90 mg/dL (0.6–5.0 mmol/L) appears. Then, press the S button to confirm the level and to go to the next step. If the test result is lower than the pre-set hypoglycemia


level, the meter will display 'Lo' as shown.

01 05-04 12:30~

▲ Caution

Step 12

Ask your healthcare professional to help you decide what your hypoglycemia level is before setting your level. Setting the Hyperglycemia (HI) Indicator



hyperglycemia indicator (possible high blood sugar). You will be alerted any time your test result is higher than the selected level. Press the ◀ or ► button until the desired hyperglycemia level between 120 and 349 mg/dL (6.7-19.4 mmol/L) appears. Press and hold the S button to confirm the hyperglycemia level and turn the meter off.



If the test result is higher than the pre-set hyperglycemia level, the meter will display $\widehat{\mathbf{w}}$ along with the measurement value. If he test result is between 241 and 600 mg/dL (13.4-33.3 mmol/L), 'KET' will blink thress times on the screen.

▲ Caution

Ask your healthcare professional to help you decide what your hyperglycemia level is before setting your level.

Note

- If your blood glucose is measured to be more than 240 mg/dL (13.3 mmol/L), you may also have ketones.
- The CareSens H Beat BT meter is not intended to detect ketones.
- If you see 'KET' displayed on your meter, ask your healthcare professional when and how you should test for ketones.

🔎 Note

If the test result is within the selected normal blood glucose range, the smile symbol will be displayed as

05-04 12:30,

10 Pairing your Meter with Smartphone

Step 1

Smartphone: If you want to pair (connect) your meter to your smartphone, launch the SmartLog mobile app and find the Accessories menu on your smartphone. Step 2

Smartphone: Select the CareSens H Beat BT model from the meters list, and then tap Bluetooth Register > Next > Search. Step 3



ا آ

Meter: With the meter turned off, press and hold the ▶ button for three seconds to enter pairing mode. The Bluetooth symbol and 'bT' will appear and 'YES' will blink at the bottom of the screen.

Step 4 Meter: Press the S button to select 'YES'. The Bluetooth symbol will blink on the screen.

Smartphone: Your smartphone will search for the meter. Tap your meter's ID (CareSens XXXX) when it appears on the Search window.

number



Step 5

Make sure the four digits of the meter device ID matches the last four digits of the meter serial number, which can be found on the meter label.



Smartphone/Meter: The Bluetooth Pairing Request window will pop up on the SmartLog mobile app screen, and a six digit PIN number will appear on the meter screen.


12.34 56 O Note

Make sure the PIN on your smartphone matches the PIN on your meter screen. Step 7

Smartphone: Tap Pair on the Bluetooth Pairing Request window.



Meter: When your meter and smartphone are successfully paired, the meter will display 'SYNC'. The saved test results will be automatically transferred to SmartLog.



End O Note

If the data transfer fails, the meter will display 'FAIL' and it will be automatically turned off after three seconds. Then, repeat steps 1 to 7.

Note

Some smartphones, especially those that are not tested or approved by i-SENS, may be incompatible with your meter. Visit www.i-sens.com/ smartlog for more information about supported smartphones. You can also scan the QR code on the back cover of this user manual.

O Note


- If your meter and smartphone are paired, the meter communicates with the smartphone, automatically synchronizing the current date and time. To set the date and time manually, press and hold the S button to enter the SET
- mode, and continue to press the **S** button to go to Adjusting the Date and Time.


11 Setting the Strip Expiration Date Indicator

Step 1 Entering the Expiration Date Setting Press and hold the ◀ and ► buttons at the same time for

3 seconds to enter the expiration date settings. After all segments flash across the screen, 'EP' will show up.

O Note





ranges printed on their labels. Compare the result displayed on the meter to the CareSens H Control Solution range printed on the test strip vial. Before using a new meter or a new vial of test strips, you may conduct a control solution test following the procedure on pages 18-19.

O Note

- Use only the CareSens H Control Solution (available for purchase • Check the expiration date printed on the bottle. When you first open
- a control solution bottle, record the discard date (date opened plus three (3) months) in the space provided on the label.
- Make sure your meter, test strips, and control solution are at room temperature before testing. Control solution tests must be done at room temperature (20–25 °C). • Before using the control solution, shake the bottle, discard the first
- Close the control solution bottle tightly and store at a temperature
- between 8-30 °C.
- You may do a control solution test:
- When you want to practice the test procedure using the control solution instead of blood,
- When using the meter for the first time,

Control Solution Testing

0504 /18³37,

15

- Whenever you open a new vial of test strips. If the meter or test strips do not function properly
- If your symptoms are inconsistent with the blood glucose test results and you feel that the meter or test strips are not working

Step 1 Inserting Test Strip

Insert a test strip into the meter's test strip

port, with the contact bars facing upwards.

Gently push the test strip into the port until the

meter beeps. Be careful not to bend the strip

while pushing it in. The been symbol will show

Step 2 Activating Control Solution Test Mode

You can flag the control solution test result by

pressing the ► button for 3 seconds. To undo

the control solution flag, press the ► button for

18

properly • If you drop or damage the meter.

3 seconds again.



Step 3 Applying Control Solution to Test Strip Shake the bottle before each test. Remove the cap and squeeze the bottle to discard the first drop. Then wipe the tip with a clean tissue or cloth. Dispense a drop of control solution onto a clean non-absorbent surface. It helps to squeeze a drop onto the top of the cap as shown. After the been symbol appears on the display, apply the solution to the edge of the test strip until the meter beeps. Make sure the confirmation window fills completely.

O Note

The meter may switch off if the control solution sample is not applied within 2 minutes of the **barren** symbol appearing on the screen. If the meter turns off, remove the strip, reinsert, and start from step 1.



Step 4 Waiting for the Result The display segments will rotate clockwise and a test result will appear after the meter counts down from 5 to 1. When flagged, the result is stored in the meter's memory but it is not included in 0504 12:37., 0504 12:37., the averages.

Step 5 Comparing the Result Compare the result displayed on the meter to the range printed on the test Control Solution Hange Control A: XX-XX mg/dL (X.X-XX mmol/L) Control B: XX-XX mg/dL (X.X-XX mmol/L) the range. strip vial. The result should fall within

68

Control Solution Range

▲ Caution

The range printed on the test strip vial is for the CareSens H Control Solution only. It has nothing to do with your blood glucose level.

Note

The CareSens H Control Solution can be purchased separately. Please contact your authorised i-SENS sales representative.

Comparing the Control Solution Test Results


The test result of each control solution should be within the range printed on the label of the test strip vial. Repeat the control solution test if the test result falls outside of the range. Out of range results

may occur in following situations.	
Situations	Do This
 When the control solution bottle was not shaken well, When the meter, test strip, or the control solution were exposed to high or low temperatures, When the first drop of the control solution was not discarded or the tip of the bottle was not wiped clean, When the meter is not functioning properly. 	Repeat the control solution test by referring to the Note on page 18.
 When the control solution is past the expiration date printed on the bottle, When the control solution is past its discard date (the date the bottle was opened plus three (3) months), When the control solution is contaminated. 	Discard the used control solution and repeat the test using a new bottle of control solution.

• when the control solution is contaminated. or control solution. If results continue to fall outside the range printed on the test strip vial, the CareSens H test strip and CareSens H Beat BT meter may not be working properly. Do not use your system and contact i-SENS sales representative.

13 Using the Lancing Device

You will need a lancing device in order to collect a blood sample. You may use the lancing device that is included in the CareSens H Beat BT Blood Glucose Monitoring System or any other medically approved lancing device.



• The lancing device is for use by a single user only and should not be shared with anyone.

• Use a soft cloth or tissue to wipe the lancing device. If necessary, a small amount of alcohol on a soft cloth or tissue may be used.

▲ Caution

To avoid infection when drawing a sample, do not use a lancet more

than once, and: • Do not use a lancet that has been used by others.

Always use a new sterile lancet.

Keep the lancing device clean.

Note


Repeated puncturing at the same sample site may cause pain or skin calluses (thick hard skin). Choose a different site each time you test.

Preparing the Lancing Device Step 1



Wash hands and sample site with soap and warm water. Rinse and dry thoroughly.





Firmly insert a new lancet into the lancet holder. Hold the lancet firmly. Gently twist to pull off protective disk. Save disk to recap

Step 4 Turn the adjustable tip until it is aligned with the load confirmation window and release button as shown.

lancet after use. Replace lancing device tip.



Note

0 = a shallow puncture for softer skin 5 = a deeper puncture for thick or calloused skin

Step 5



To cock the lancing device, hold the body of lancing device in one hand and pull the sliding barrel with the other hand. The device is loaded when you feel a click and the load

confirmation window turns red.

The lancing device has six puncture depth

settings (0 for a shallow puncture, 5 for a

desired number aligns with the arrow.

deeper puncture). Choose a depth by rotating

the top portion of the adjustable tip until the

🔎 Note

-0-7

The skin depth to get blood samples will vary for various people at different sample sites. The lancing device's adjustable tip allows the best depth of skin penetration to get an adequate sample size.

Preparing the Meter and Test Strip Step 7

Insert a test strip with the contact bars facing upwards into the meter's test strip port. Push the strip in gently until the meter beeps. Be careful not to bend the test strip. The **equal** symbol will appear on the screen.



Step 8





Obtain a blood sample using the lancing device. Place the device against the pad of the finger. The best puncture sites are on the middle or ring fingers. Press the release button. Remove the device from the finger. Wait a few seconds for a blood drop to form. A minimum volume of 0.5 microliter is needed to fill the confirmation window (actual size of 0.5 µL: ●).

Step 9

After the best symbol appears on the screen, apply the blood sample to the end of the test strip till the meter beeps. If the confirmation window is not filled in time because of abnormal viscosity (thickness and stickiness) or insufficient volume, the Er4 message may appear. It is recommended to place the test strip vertically into the blood sample site as shown below.



▲ Caution

Do not allow any foreign substances, such as dirt, blood, or water, enter into the meter. The meter may be damaged or may malfunction. Follow the warning information provided below to prevent possible damage to the meter. • Do not apply the blood sample directly to the test strip port.

• Do not apply the blood sample to the test strip while holding the meter in a way that the tip of the test strip faces upwards. The blood sample may run down the surface of the test strip and flow into the test strip port. Do not store your meter in unsanitary or contaminated sites.

🔎 Note

The meter may switch off if the blood sample is not applied within 2 minutes of the bar symbol appearing on the screen. If the meter turns off, remove the strip and reinsert it, and start from Step 2. 23

Step 10

Apply the blood sample to the end of the test strip until you hear a 'beep'. At this time, the display segments will rotate clockwise while the blood is aoina in.

The test result will appear after the meter counts down from 5 to 1. The result will be automatically stored in the meter's memory. If the test strip is removed after the test result is displayed, the meter will automatically switch off after 3 seconds. Discard used test strips safely in disposable containers. If the Bluetooth feature is activated, the meter will send the test result to the connected smartphone.



data exists).

Note

To transmit glucose data using the Bluetooth feature, • The Bluetooth feature on the meter must be turned on, • The meter and a smartphone must be paired, The SmartLog mobile app must be launched. The meter will transmit data in the following cases, • When the strip is ejected after measuring, • When the meter is turned on (only when untransmitted

Step 11

21

05-04 12:37...


You can attach a flag to a result to indicate particular situations while the strip is still in the meter. When the result is displayed right after a test, press the ◀ or ► button to select a pre-meal flag (🍎), a postmeal flag(), a fasting flag(). When you remove the test strip while the desired flag is blinking, the test result is stored with the flag. If you do not want to add any flags on the test result, remove the strip after the test result is displayed.



24

05-04 12:37... 0504 12:37... 05-04 I2:31~ 05:04 12:37... Post-meal flag Fasting flag No flac Pre-meal flag Discarding Used Lancets Step 1 Unscrew the lancing device tip.






Step 2



The lancet is for single use only. Never share or reuse a lancet. Always

dispose of lancets properly.

14 HI and Lo Messages



Lo Message 'Lo' appears when a test result is less than 10 mg/dL (0.6 mmol/L) and indicates severe hypoglycemia (very low glucose levels).

If 'Lo' is displayed again upon retesting, please contact US-DH 12:37... your healthcare professional immediately.

O Note

22

sample

Please contact your authorised i-SENS sales representative if such messages are displayed even though you do not have hyperglycemia or hypoglycemia.

15 Target Blood Glucose Ranges

Reminders Time of day	Your target ranges from your healthcare professional		
Before breakfast			Note
Before lunch or dinner			control
1 hour after meals		WILLI	U Syrric
2 hours after meals			
Between 2 a.m. and 4 a.m.		18	3 Sett


Expected Values

Normal blood glucose levels for an adult without diabetes are below 100 mg/dL (5.5 mmol/L) before meals and fasting* and are less than 140 mg/dL (7.8 mmol/L) two hours after meals. * Fasting is defined as no caloric intake for at least eight hours.

Reference

American Diabetes Association (Standards of Medical Care in Diabetes - 2021. Diabetes Care), January 2021, vol. 44 (Supplement 1): S15-S33.

16 Transferring Test Results Using Cable




05-04 12:30... www.i-sens.com. O Note

t results can also be transferred wirelessly using Bluetooth. To pair the meter and your smartphone, see page 15.


17 Meter Memory

The CareSens H Beat BT meter can save up to 1,000 glucose test results with time and date. If the memory is full, the oldest test result will be deleted and the latest test result will be stored. The meter calculates and displays the averages of total test results, pre-meal (\bigstar) test results, post-meal test (Υ), and fasting test results (() from the last 1, 7, 14, 30 and 90 days.

Viewing Averages Stored in Memory



The number of taste of tests within the displayed at the bottom of the screen followed current day by the 1 day average value and the number of the test results saved within the current day. 27




Step 6 Press the S button to turn off the meter.

25

O Note in the averages.

of tests within the current day





(On) 05: 10-

O Note

0.4

797

(OFF) OS: 10~

0

ঠ12:00∞

0

Step 1

Step 2



7, 14, 30 and 90-day average Press the < button to view 1. values and the number of tests performed for the last test



period. Repeatedly press the ◀ button

with the **š** symbol for the last test period.



7, 14, 30 and 90-day average values and the number of tests performed post-meals with the Y symbol for the last test Press the
< button to view 1. x symbol for the last test

Press the ◀ button to view 1, 7, 14, 30 and 90-day average values and the number of tests performed during fasting with the (?) symbol for the last test values and the number of tests the (?) symbol for the last test 9077725 5n period.

28

Use the ▶ button to scroll back through the averages seen previously.

The control solution test results saved with the **b** symbol are not included

Viewing Test Results Stored in Memory

Step 1

The number Press the \triangleleft or \blacktriangleright or \mathbf{S} button to turn the meter on. The current date and time will be displayed on the bottom of the screen followed by the 1 day average value and the number of the test results saved within the current day.



Step 2 Use the ► button to scroll through the test results, starting from the most recent and ending with the oldest. Press the < button to return to 08-03 ... 12:03-

the result seen previously. After checking the stored test results hold the **S** button to turn off the meter.

The control solution test results saved with a symbol will be displayed with a symbol when you review the stored test results.

18 Setting the Alarm Function

Four types of alarms can be set in the CareSens H Beat BT meter: one post-meal alert (PP2 alert) and three time set alarms (alarm 1-3). The PP2 alert goes off 2 hours after setting the alert. The alarms ring for 15 seconds and can be silenced by pressing the ◀ or ► or S button or by inserting a test strip. 29

Setting the Post-meal Alert (PP2 alert)



Step 1 Turning the PP2 alert On Without inserting a test surp, press and hold the button for 3 seconds to set the post-meal alert. 'PP2', bell (🌲) symbol and 'On' will be displayed. The screen will then automatically change to the memory recall mode. At this time, bell (🌲) symbol, indicating that the PP2 alert has been set. will be displayed on the screen.

The PP2 alert will automatically turn off if the meter's time setting is adjusted to more than two hours before or just past the currently activated PP2 alert time.



Step 2 Turning the PP2 alert OFF To turn off the PP2 alert, press and hold the ◄ button for 3 seconds. PP2', bell (♠) symbol and 'OFF' will appear on the screen. Then the screen will change automatically to the memory recall mode without bell - A disappears (A) symbol displayed.

Setting the Time Alarms (alarm 1-3)

Without inserting a test strip, press the < and S buttons simultaneously for 3 seconds to enter the time alarm setting, 'alarm 1' will be displayed while 'OFF' blinks on the screen.

On pressing the ▶ button, 'alarm 1' is set and 'On' is displayed on the screen. Press the ► button again to cancel 'alarm 1'. 'OFF' will blink on the screen.


30

Step 3

0

Press the < button to adjust the time of 'alarm 1'. A number representing the hour will blink on the screen. Press the button to set the hour.

0n di:00,, Step 4



Press the S button to finish and to go to 'alarm 2' settng Repeat steps 2 to 4 to set the remaining time alarms

(alarm 2–3).

0. 240 - 12:00 - 12:00 - 12:00 - 12:00 - 12:00 - 12:00 - 12:00 - 12:00 - 12:00 - 12:00 - 12:00 - 12:00 - 12:00 - 1 Step 6

Press the S button for 3 seconds to finish and turn the meter off.

19 Understanding Error Messages

Message	What It Means	What To Do
Erl	A used test strip was inserted.	Repeat the test with a new test strip.
		E

8-2	The blood or control solution sample was applied before the beam symbol appeared.	Repeat the test with a new test strip and wait until the symbol appears before applying the blood or control solution sample.
Er 3	The temperature during the test was above or below the operating range.	Move to an area where the temperature is within the operating range $(5-45^\circ\text{C})$ and repeat the test after the meter and test strips have reached a temperature within the operating range.
ErY	The blood sample has abnormally high viscosity or insufficient volume.	Repeat the test after inserting a new test strip.
ErS	This error message may appear when the wrong blood glucose test strip is used instead of CareSens H blood glucose test strip.	Repeat the test with a CareSens H test strip.
Erb	There is a problem with the meter.	Do not use the meter. Contact your authorised i-SENS sales representative.

Ern	There is a problem with Bluetooth communication	Contact your authorised i-SENS sales representative.
Er8	An electronic error occurred during the test.	Repeat the test with a new test strip. If the error message persists, contact your authorised i-SENS sales representative.

O Note

If the error messages persist, contact your authorised i-SENS sales representative.

20 General Troubleshooting

Problem	Troubleshooting
The display is blank even after inserting a test strip.	 Check whether the test strip is inserted with the contact bars facing up. Check if the strip has been inserted completely into the test strip port. Check if the appropriate test strip was used. Check whether the the battery is inserted with the + side facing up. Replace the battery.
The test does not start even after applying the blood sample on the strip.	Check if the confirmation window is filled completely.Repeat the test after inserting a new test strip.
The test result doesn't match the way you feel.	 Repeat the test after inserting a new test strip. Check the expiration date of the test strip. Perform control solution test.

O Note

If the problem is not resolved, please contact your authorised i-SENS sales representative.

21 Performance Characteristics

System Accuracy and Measurement Precision The performance of CareSens H Beat BT Blood Glucose Monitoring System has been evaluated in laboratory and in clinical tests. Accuracy: The accuracy of the CareSens H Beat BT Blood Glucose Monitoring System (Model GM01YAB) was assessed by comparing blood glucose results obtained by patients with those obtained using a YSI Model 2300 Glucose Analyzer, a laboratory instrument. The following results were obtained by diabetic patients at clinic centers.

0		, , , , , , , , , , , , , , , , , , , ,		
Slope		0.980		
Y-intercept		3.682 mg/dL (0.2	0 mmol/L)	
Correlation coefficient (r	.)	0.9952		
Number of Tests		100		
Range tested		31.6-458.0 mg/d	L (1.8–25.4 mmol/L)	
System accuracy results for	or gluco	ose concentration <	: 100 mg/dL (5.55 mmol/L)	
Within ±5 mg/dL (Within ±0.28 mmol/L)		hin ±10 mg/dL in ±0.56 mmol/L)	Within ±15 mg/dL (Within ±0.83 mmol/L)	
130/180 (72.2 %)	173/180 (96.1 %)		177/180 (98.3 %)	
System accuracy results for glucose concentration ≥ 100 mg/dL (5.55 mmol/				
Within ±5 %	V	Vithin ±10 %	Within ±15 %	
301/420 (71.7 %)	398	3/420 (94.8 %)	415/420 (98.8 %)	
System accuracy results for glucose concentrations between 31.6 mg/dL (1.8 mmol/L) and 458 mg/dL (25.4 mmol/L)				
Within ±15 mg/dL (Within ±0.83 mmol/L) and Within ±15 %				
592/600 (98.7 %)				
D	1		anne and the labele eventages	

Precision: The precision studies were performed in a laboratory using CareSens H Beat BT Blood Glucose Monitoring Systems.

Within Run Pre	ecision		
	38 mg/dL (2.1 mmol/L)	SD = 2.0 mg/dL (0.1 mmol/L)	
	83 mg/dL (4.6 mmol/L)	SD = 3.6 mg/dL (0.2 mmol/L)	
Blood average	130 mg/dL (7.2 mmol/L)	CV = 3.2 %	
	192 mg/dL (10.7 mmol/L)	CV = 2.8 %	
	312 mg/dL (17.3 mmol/L)	CV = 2.5 %	
Between Run Precision			
Control	35 mg/dL (2.0 mmol/L)	SD = 1.5 mg/dL (0.1 mmol/L)	
Control	120 mg/dL (6.7 mmol/L)	CV = 4.0 %	
average	347 mg/dL (19.3 mmol/L)	CV = 4.1 %	

Influence Quantities

Packed Cell Volume (Hematocrit)

Packed cell volume evaluation was conducted in various hematocrit levels. The range of hematocrit levels within the acceptance criteria is 15–70 %. Interferences

The effect of various interfering substances was evaluated in whole blood samples. The presence of the following substances within the given concentrations does not affect blood glucose measurements. Higher concentrations of the substances shown below may cause inaccurate blood glucose results.

No.	Interferent	Concentration	No.	Interferent	Concentration
1	Acetaminophen (paracetamol)	20 mg/dL (1.32 mmol/L)	13	Ibuprofen	50 mg/dL (2.42 mmol/L)
2	Ascorbic acid	3 mg/dL (0.17 mmol/L)	14	lcodextrin	1094 mg/dL
3	Bilirubin	20 mg/dL (0.34 mmol/L)	15	L-Dopa (L-3,4- dihydroxyphenylalanine)	5 mg/dL (0.25 mmol/L)
4	Cholesterol	500 mg/dL (12.93 mmol/L)	16	Maltose	1000 mg/dL (29.21 mmol/L)
5	Creatinine	30 mg/dL (2.65 mmol/L)	17	Methyl-DOPA	1.5 mg/dL (0.07 mmol/L)
6	Dopamine	13 mg/dL (0.85 mmol/L)	18	Pralidoxime Iodide (PAM)	25 mg/dL (0.95 mmol/L)
7	EDTA	180 mg/dL (6.16 mmol/L)	19	Sodium Salicylate	70 mg/dL (5.07 mmol/L)
8	Galactose	60 mg/dL (3.33 mmol/L)	20	Tolbutamide	100 mg/dL (3.70 mmol/L)
9	Gentisic acid	50 mg/dL (3.24 mmol/L)	21	Tolazamide	100 mg/dL (3.21 mmol/L)
10	Glutathione (Red)	92 mg/dL (2.99 mmol/L)	22	Triglycerides	3300 mg/dL (37.26 mmol/L)
11	Hemoglobin	500 mg/dL (0.31 mmol/L)	23	Uric acid	25 mg/dL (1.49 mmol/L)
12	Heparin	8000 U/dL	24	Xylose	7.2 mg/dL (0.48 mmol/L)

Compounds of xylose \geq 7.2 mg/dL (0.48 mmol/L) at glucose concentrations of 50-100 mg/dL (2.78 mmol/L-5.56 mmol/L) may cause overestimation of blood glucose results.

User Performance Evaluation

32

A study evaluating glucose values from fingertip capillary blood samples obtained by 100 lay persons showed the following results:

100 % within ±15 mg/dL (±0.83 mmol/L) of the medical laboratory values at glucose concentrations below 100 mg/dL (5.55 mmol/L), and 100 % within ±15 % of the medical laboratory values at glucose concentrations at or above 100 mg/dL (5.55 mmol/L).

22 Warranty Information

Manufacturer's Warranty

i-SENS, Inc. warrants that the CareSens H Beat BT meter shall be free of defects in material and workman a period of five (5) years. The meter must have been subjected to normal use. The warranty does not cover improper handling, tampering, use, or service of the meter. Any claim must be made within the warranty period.

i-SENS will, at its discretion, repair or replace a defective meter or meter part that is covered by this warranty. As a matter of warranty policy, i-SENS will not reimburse the consumer's purchase price.

Obtaining Warranty Service

To obtain warranty service, you must return the defective meter or meter part along with proof of purchase to your nearest i-SENS sales or customer service representative.

i-SENS, Inc. 43, Banpo-daero 28-gil Seocho-gu, Seoul 06646, Korea www.i-sens.com	EC REP Medical Technology Promedt Consulting GmbH Altenhofstrasse 80 66386 St. Ingbert, Germany
---	---

C € 0123 **IVD**

©2023 i-SENS, Inc. All Rights Reserved. 36